

! TFP Ø 18 mm  
no longer available!

✓ Successor: TSMF  
New, modular and better!

All advantages at [anderson-negele.com](http://anderson-negele.com)



HYGIENIC BY DESIGN

ANDERSON-NEGELE

Product Information TFP-641, -642, -661, -681, -841, -842, -861, -881

PHARMA

# Temperature Sensor PHARMadapt EPA



## Application / Specified Usage

- Developed for applications in pharmaceutical industry and biotechnology
- Temperature measurement especially at small pipe diameters
- In connection with build-in system **PHARMadapt EPA** suitable for pipes DN10...100

## Application Examples

- Process monitoring
- Monitoring of CIP-/ SIP-cleaning

## Hygienic Design / Process Connection

- Hygienic and easy sterilizable installation by using Negele build-in system **PHARMadapt EPA**
- CIP-/ SIP-cleaning up to 140 °C
- All product contacting materials compliant to FDA
- Sensor completely made of stainless steel
- Sealing ring according to USP Class VI
- Conforming to 3-A Sanitary Standard for DIN 11866 series A with DN ≥ 25, DIN 11866 series B with DN ≥ 20, DIN 11866 series C with DN ≥ 1"

## Features

- For small pipe diameters from DN10
- Easy demounting for cleaning and calibration by clamp system
- Short reaction time, very compact measure point with leakage control
- Integrated transmitter available
- Light weight sensor head, non-sensitive to vibrations
- Hygienic lid design
- Electrical connection via M12 plug
- Material 1.4435, material certificate 3.1 in scope of delivery (for all product contacting parts)
- Quick and easy to install with an orbital welding machine


## Options / Accessories

- 2x Pt100, optional (not retrofittable)
- 2x Pt100 with two transmitters (not retrofittable)
- Programmable transmitter MPU-4 and MPU-M with 4...20 mA output, 2-wire
- Integrated transmitter for HART-protocol
- Programming adapter MPU-P 9701
- Integrated display MPU-LCD in connecting head
- Pt100-chip with other classes of accuracy (1/3 B, 1/10 B)
- Preassembled cable for M12 plug
- Fixed cable in other length or material available
- Calibration certificate (just available with order placement)
- Customer specific label and TAG-number plate (stainless steel)

## Authorisations



## Temperature sensor TFP-641 with build-in system EPA-18



## Temperature sensor TFP-661 with transmitter MPU-M and TAG-number plate



<b>Temperature sensor</b>		
<b>Process connection</b>	gap free	with clamp-ring SRC-05 resp. SRC-10
<b>Insertion length EL</b>	TFP-6xx TFP-8xx	10 mm, 25 mm, 50 mm, 100 mm 20 mm, 50 mm
<b>Materials</b>	connection head protection tube sealing ring	stainless steel 1.4301 (AISI 304) stainless steel 1.4435 (AISI 316L) EPDM, USP Class VI, FDA 21 CFR 177.2600
<b>Temperature ranges</b>	ambient sensor tip	-50...+80 °C -50...+250 °C
<b>Operating pressure</b>		10 bar max.
<b>Sensing resistor</b>	acc. to DIN EN 60751	Pt100
<b>Electrical connection</b>	TFP-64x, TFP-84x  TFP-661, TFP-861 TFP-681, TFP-881	cable gland M16x1.5 (PG) or M12 plug 1.4301 (AISI 304), 4-pin M12 plug 1.4301 (AISI 304) fixed cable (PTFE, 4x 0.14 mm <sup>2</sup> ), standard: 2.5 m
<b>Protection type</b>		IP 69 K (with electrical connection M12 plug)

<b>Transmitter MPU-4, MPU-H, MPU-M</b>		
<b>Temperature ranges</b>	ambient storage	-40...+85 °C -55...+90 °C
<b>Measuring ranges</b>	MPU-4, MPU-H, MPU-M	standard: -10...40 °C, 0...50 / 100 / 150 / 200 °C special ranges free programmable
<b>Accuracy</b>	input	< ±0.25 °C
<b>Temperature drift</b>	zero, span	< 0.01 % / K
<b>Supply</b>	MPU-4, MPU-H, MPU-M accuracy	8...35 V DC 0.01 % / V (reference: 12 V DC)
<b>Output</b>	signal accuracy burden	analog 4...20 mA < ±0.1 % of measurement range < 600 Ω (at U <sub>B</sub> = 24 V)
<b>Humidity</b>	without condensation	0...98 %

**Accuracy classes of temperature sensors | Tolerances for Pt100 acc. to DIN EN 60751**

Pt100	A	1/3 B	1/10 B
0 °C / 100 Ω	±0.15 K / ±0.06 Ω	±0.10 K / ±0.04 Ω	±0.03 K / ±0.01 Ω
100 °C / 138.5 Ω	±0.35 K / ±0.13 Ω	±0.27 K / ±0.10 Ω	±0.08 K / ±0.03 Ω


<b>Table reaction time</b>	<b>EPA-8, EPA-18</b>
t <sub>50</sub>	4.4 s
t <sub>90</sub>	13.1 s

**Reaction time**



The mentioned times were measured by emersing a temperature sensor from room temperature into boiling water.

Electrical connection without transmitter

With 1x M12 plug




Configuration 1st M12 plug



Electrical connection with transmitter


With M12 plug




Configuration M12 plug

- 1: + supply
- 2: - supply 4...20 mA
- 3: not connected
- 4: not connected


With 2x M12 plug




Configuration 2nd M12 plug




With cable gland



With cable gland




Configuration strip terminal



- clamps for 1st Pt100
- clamps for 2nd Pt100 (at version 2x Pt100)

Electrical connection with two transmitter (TFP-642, -842)

With 1x M12 plug (sensor 1 + sensor 2)



Configuration M12 plug


- 1: + supply (sensor 1)
- 2: - supply 4...20 mA (sensor 1)
- 3: - supply 4...20 mA (sensor 2)
- 4: + supply (sensor 2)

With fixed cable




Fixed cable connection with 1x Pt100

wh ye bn gn standard  
rd rd wh wh PTFE




With 2x M12 plug (sensor 1)



Configuration M12 plug

- 1: + supply (sensor 1)
- 2: - supply 4...20 mA (sensor 1)
- 3: not connected
- 4: not connected


With 2x M12 plug (sensor 2)




Configuration M12 plug

- 1: + supply (sensor 2)
- 2: - supply 4...20 mA (sensor 2)
- 3: not connected
- 4: not connected


TFP-641 | TFP-641.2 with cable gland




TFP-642 with M12 plug




TFP-661




TFP-661 / ... / MPU-M




TFP-681




TFP-841 | TFP-841.2 with cable gland




TFP-842 with M12 plug




TFP-861



TFP-861 / ... / MPU-M



TFP-881



**Mechanical Connection / Installation**



- The sensors are only for use with pharmaceutical build-in system **PHARMadapt EPA**.

**Mounting Advice**



- The suitable insertion length depends on the pipe diameter of the measuring point. See dimension tables in product information **PHARMadapt EPA**.

**Transport / Storage**



- No outdoor storage
- Dry and dust free
- Not exposed to corrosive media
- Protected against solar radiation
- Avoiding mechanical shock and vibration
- Storage temperature -55...+90 °C
- Relative humidity maximum 98 %

**Conventional Usage**



- Not suitable for applications in explosive areas.
- Not suitable for applications in security-relevant equipments (SIL).

**Reshipment**



- Sensors shall be clean and must not be contaminated with dangerous media!
- Use suitable transport packaging only to avoid damage of the equipment!

**Cleaning / Maintenance**



- In case of using pressure washers, don't point nozzle directly to electrical connections!

**Standards and Guidelines**



- You have to comply with applicable regulations and directives.

**Disposal**



- Electrical devices should not be disposed of with household trash. They must be recycled in accordance with national laws and regulations.
- Take the device directly to a specialized recycling company and do not use municipal collection points.

**Note on CE**



- Applicable directives: Electromagnetic Compatibility Directive 2014/30/EU
- Compliance with the applicable EU directives is identified by the CE label on the product.
- The operating company is responsible for complying with the guidelines applicable to the entire installation.

**Conditions for a measuring point according to 3-A Sanitary Standard 74-06**



- The sensors TFP-641, -642, -661, -681, -841, -842, -861, -881 conforming to the 3-A Sanitary Standard.
- The sensors are designed for CIP-/ SIP-cleaning. Maximum 140 °C / 120 minutes.
- Only with the build-in system **PHARMadapt EPA** allowed.
- Mounting position, self draining and the position of the leakage hole must be in accordance to current 3-A Sanitary Standard.

**Accessories**

**PVC-cable with M12-connection made of 1.4305, IP 69 K, unshielded**  
**M12-PVC/4**      PVC-cable 4-pin, length 5 m, 10 m, 25 m

**PVC-cable with M12-connection, brass nickel-plated, IP 67, shielded**  
**M12-PVC/4G**      PVC-cable 4-pin, length 5 m, 10 m, 25 m

**Programming adapter**  
**MPU-P 9701**      Programming adapter for MPU-4, MPU-H and MPU-M

**Sealing ring**  
**DRE-5**      Sealing ring for EPA-8, Ø 5 x 1.5 mm, material EPDM (FDA compliant, USP class VI)  
**DRE-15**      Sealing ring for EPA-18, Ø 15 x 1.5 mm, material EPDM (FDA compliant, USP class VI)

## Order code for version with 1x Pt100

<b>TFP-641</b>	(for PHARMadapt EPA-8, connection head Ø 49 mm, non-sensitive to vibrations)
<b>TFP-661</b>	(for PHARMadapt EPA-8, connection head Ø 18 mm, electrical connection via M12 plug)
<b>TFP-681</b>	(for PHARMadapt EPA-8, connection head Ø 18 mm, electrical connection via 2.5 m PTFE-cable)
<b>TFP-841</b>	(for PHARMadapt EPA-18, connection head Ø 49 mm, non-sensitive to vibrations)
<b>TFP-861</b>	(for PHARMadapt EPA-18, connection head Ø 18 mm, electrical connection via M12 plug)
<b>TFP-881</b>	(for PHARMadapt EPA-18, connection head Ø 18 mm, electrical connection via 2.5 m PTFE-cable)

## Sensor Length EL for TFP-6xx in mm

<b>010</b>	(length 10 mm)
<b>025</b>	(length 25 mm)
<b>050</b>	(length 50 mm)
<b>100</b>	(length 100 mm)

## Sensor Length EL for TFP-8xx in mm

<b>020</b>	(length 20 mm)
<b>050</b>	(length 50 mm)

## Accuracy Class Pt100

**A**  
**1/3B**  
**1/10B**

## Electrical Connection for TFP-641 and TFP-841

<b>PG</b>	(cable gland M16x1.5)
<b>M12</b>	(M12 plug, standard with MPU-LCD)

## Transmitter

**X** (without)

## Transmitter for TFP-641 and TFP-841

<b>MPU-4</b>	(programmable)
<b>MPU-H</b>	(HART-protocol)
<b>MPU-LCD</b>	(with display)

## Transmitter for TFP-661 and TFP-861

**MPU-M** (programmable)

## Measurement Range(only for types with transmitter; not selectable at MPU-LCD)

<b>-10...40</b>	(measuring range -10...+40 °C)
<b>0...50</b>	(measuring range 0...+50 °C)
<b>0...100</b>	(measuring range 0...+100 °C)
<b>0...150</b>	(measuring range 0...+150 °C)
<b>0...200</b>	(measuring range 0...+200 °C)
<b>xx...yy</b>	(special range)



**TFP Ø 18 mm**  
no longer available!



**Successor: TSMF**  
New, modular and better!

All advantages at [anderson-negele.com](http://anderson-negele.com)

**TFP-641 / 025 / A / M12 / MPU-4 / 0...100**

## Note

The clamp ring is not included in scope of delivery and must be ordered separately.

<b>SRC-05</b>	Clamp-tension ring for EPA-8, material 1.4301
<b>SRC-10</b>	Clamp-tension ring for EPA-18, material 1.4301

## Build-In systems

Suitable build-in systems for temperature sensors TFP-641, -642, -661, -681, -841, -842, -861, -881 you will find in product information **Process Connection PHARMadapt EPA**.

Order code for version with 2x Pt100

- TFP-641.2** (for PHARMadapt EPA-8, connection head Ø 49 mm, 2x Pt100, non-sensitive to vibrations, no transmitter possible!)
- TFP-642** (for PHARMadapt EPA-8, like TFP-641.2, but higher connection head and prepared for 2x transmitter)
- TFP-841.2** (for PHARMadapt EPA-18, connection head Ø 49 mm, 2x Pt100, non-sensitive to vibrations, no transmitter possible!)
- TFP-842** (for PHARMadapt EPA-18, like TFP-841.2, but higher connection head and prepared for 2x transmitter)

**Sensor Length EL for TFP-6xx in mm**

- 010** (length 10 mm)
- 025** (length 25 mm)
- 050** (length 50 mm)

**Sensor Length EL for TFP-8xx in mm**

- 020** (length 20 mm)
- 050** (length 50 mm)

**Accuracy Class Pt100**

- A**
- 1/3B**
- 1/10B**

**Electrical Connection for TFP-641.2 and TFP-841.2**

- PG** (cable gland M16x1.5)
- 2xPG** (2x cable gland M16x1.5)
- 2xM12** (2x M12 plug)

**Electrical Connection for TFP-642 and TFP-842**

- M12** (M12 plug)
- 2xM12** (2x M12 plug)

**Continue if TFP-642 or TFP-842 is selected!  
No further options for TFP-641.2 or TFP-841.2!**

**1. Transmitter**

**MPU-4** (programmable)

**Measurement Range 1. MPU**



- 10...40** (measurement range -10...+40 °C)
- 0...50** (measurement range 0...+50 °C)
- 0...100** (measurement range 0...+100 °C)
- 0...150** (measurement range 0...+150 °C)
- 0...200** (measurement range 0...+200 °C)
- xx...yy** (special range)

**2. Transmitter**

**MPU-4** (programmable)

**Measurement Range 2. MPU**

- 10...40** (range -10...+40 °C)
- 0...50** (range 0...+50 °C)
- 0...100** (range 0...+100 °C)
- 0...150** (range 0...+150 °C)
- 0...200** (range 0...+200 °C)
- xx...yy** (special range)

 <p><b>TFP Ø 18 mm</b> no longer available!</p>	 <p><b>Successor: TSMF</b> New, modular and better!</p>
<p>All advantages at <a href="http://anderson-negele.com">anderson-negele.com</a></p>	

**TFP-642 / 025 / A / M12 / MPU-4 / 0...50 / MPU-4 / 0...50**