HYGIENIC BY DESIGN

Product Information TFP-641, -642, -661, -681, -841, -842, -861, -881

Temperature Sensor PHARMadapt EPA

Application / Specified Usage

- · Developed for applications in pharmaceutical industy and biotechnology
- · Temperature measurement especially at small pipe diameters
- In connection with build-in system PHARMadapt EPA suitable for pipes DN10...100

Application Examples

- · Process monitoring
- · Monitoring of CIP-/ SIP-cleaning

Hygienic Design / Process Connection

- Hygienic and easy sterilizable installation by using Negele build-in system PHARMadapt EPA
- · CIP-/ SIP-cleaning up to 140 °C
- · All product contacting materials compliant to FDA
- · Sensor completely made of stainless steel
- · Sealing ring according to USP Class VI
- Conforming to 3-A Sanitary Standard for DIN 11866 series A with DN \ge 25, DIN 11866 series B with DN \ge 20, DIN 11866 series C with DN \ge 1"

Features

- · For small pipe diameters from DN10
- · Easy demounting for cleaning and calibration by clamp system
- · Short reaction time, very compact measure point with leackage control
- Integrated transmitter available
- · Light weight sensor head, non-sensitive to vibrations
- · Hygienic lid design
- · Electrical connection via M12 plug
- Material 1.4435, material certificate 3.1 in scope of delivery (for all product contacting parts)
- · Quick and easy to install with an orbital welding machine

Options / Accessories

- 2x Pt100, optional (not retrofittable)
- 2x Pt100 with two transmitters (not retrofittable)
- $\cdot\,$ Programmable transmitter MPU-4 and MPU-M with 4...20 mA output, 2-wire
- · Integrated transmitter for HART-protocol
- Programming adapter MPU-P 9701
- Integrated display MPU-LCD in connecting head
- Pt100-chip with other classes of accuracy (1/3 B, 1/10 B)
- \cdot Preassemled cable for M12 plug
- $\cdot\,$ Fixed cable in other length or material available
- · Calibration certificate (just available with order placement)
- \cdot Customer specific label and TAG-number plate (stainless steel)

Authorisations



build-in system EPA-18



Temperature sensor TFP-661 with transmitter MPU-M and TAG-number plate



PHARMA

PHARM adapt

PHARMA

2

| Temperature sensor | | |
|-----------------------|--|--|
| Process connection | gap free | with clamp-ring SRC-05 resp. SRC-10 |
| Insertion length EL | ТFР-6хх ТFР-8хх | 10 mm, 25 mm, 50 mm, 100 mm 20 mm, 50 mm |
| Materials | connection head protection tube sealing ring | stainless steel 1.4301 (AISI 304) stainless steel 1.4435 (AISI 316L) EPDM, USP Class VI, FDA 21 CFR 177.2600 |
| Temperature ranges | ambient sensor tip | -50+80 °C -50+250 °C |
| Operating pressure | | 10 bar max. |
| Sensing resistor | acc. to DIN EN 60751 | Pt100 |
| Electrical connection | TFP-64x, TFP-84x TFP-661, TFP-861 TFP-681, TFP-881 | cable gland M16x1.5 (PG) or M12 plug 1.4301 (AISI 304), 4-pin M12 plug 1.4301 (AISI 304) fixed cable (PTFE, 4x 0.14 mm²), standard: 2.5 m |
| Protection type | | IP 69 K (with electrical connection M12 plug) |

Transmitter MPU-4, MPU-H, MPU-M

| Temperature ranges | ambient storage | -40+85 °C -55+90 °C | | | | |
|--------------------|---------------------------------|--|--|--|--|--|
| Measuring ranges | MPU-4, MPU-H, MPU-M | standard: -1040 °C, 050 / 100 / 150 / 200 °C special ranges free programable | | | | |
| Accuracy | input | < ±0.25 °C | | | | |
| Temperature drift | zero, span | < 0.01 % / K | | | | |
| Supply | MPU-4, MPU-H, MPU-M accuracy | 835 V DC 0.01 % / V (reference: 12 V DC) | | | | |
| Output | signal accuracy burden | analog 420 mA < ±0.1 % of measurement range < 600 Ω (at U _B = 24 V) | | | | |
| Humidity | without condensation | 098 % | | | | |

Accuracy classes of temperature sensors | Tolerances for Pt100 acc. to DIN EN 60751

| Pt100 | Α | 1/3 B | 1/10 B |
|------------------|-------------------|-------------------|-------------------|
| 0 °C / 100 Ω | ±0.15 K / ±0.06 Ω | ±0.10 K / ±0.04 Ω | ±0.03 K / ±0.01 Ω |
| 100 °C / 138.5 Ω | ±0.35 K / ±0.13 Ω | ±0.27 K / ±0.10 Ω | ±0.08 K / ±0.03 Ω |

| Table reaction time | EPA-8, EPA-18 | Reaction time |
|---------------------|---------------|---|
| t ₅₀ | 4.4 s | The mentioned times were measured by emersing a |
| t ₉₀ | 13.1 s | temperature sensor from room temperature into boiling |
| | | water. |





+ supply (sensor 2)
- supply 4...20 mA (sensor 2)
not connected
not connected

TFP-641 | TFP-641.2 with cable gland



ø4

TFP-642 with M12 plug





TFP-841 | TFP-841.2 with cable gland



TFP-842 with M12 plug





PHARMA

Mechanical Connection / Installation

• The sensors are only for use with pharmaceutical buildin system **PHARMadapt EPA**.

Transport / Storage

- No outdoor storage
- · Dry and dust free

Reshipment

Disposal

- Not exposed to corrosive media
- · Protected against solar radiation
- · Avoiding mechanical shock and vibration

Sensors shall be clean and must not be contaminated

Electrical devices should not be disposed of with house-

hold trash. They must be recycled in accordance with

Take the device directly to a specialized recycling com-

pany and do not use municipal collection points.

Use suitable transport packaging only to avoid

Storage temperature -55...+90 °C

with dangerous media!

damage of the equipment!

national laws and regulations.

· Relative humidity maximum 98 %





• Not suitable for applications in explosive areas.

The suitable insertion length depends on the pipe

in product information PHARMadapt EPA.

diameter of the measuring point. See dimension tables

 Not suitable for applications in security-relevant equipments (SIL).

Cleaning / Maintenance

Mounting Advice

Conventional Usage

 In case of using pressure washers, dont't point nozzle directly to electrical connections!

Standards and Guidelines

You have to comply with applicable regulations and directives.

Note on CE

- · Applicable directives:
- Electromagnetic Compatibility Directive 2014/30/EU
- Compliance with the applicable EU directives is identified by the CE label on the product.
- The operating company is responsible for complying with the guidelines applicable to the entire installation.
- Conditions for a measuring point according to 3-A Sanitary Standard 74-06



- \cdot The sensors TFP-641, -642, -661, -681, -841, -842, -861, -881 conforming to the 3-A Sanitary Standard.
- The sensors are designed for CIP-/ SIP-cleaning. Maximum 140 °C / 120 minutes.
- \cdot Only with the build-in system <code>PHARMadapt EPA</code> allowed.
- Mounting position, self draining and the position of the leackage hole must be in accordance to current 3-A Sanitary Standard.

Accessories

| PVC-cable with M1 | 2-connection made of 1.4305, IP 69 K, unshielded |
|-------------------|---|
| M12-PVC/4 | PVC-cable 4-pin, length 5 m, 10 m, 25 m |
| PVC-cable with M1 | 2-connection, brass nickel-plated, IP 67, shielded |
| M12-PVC/4G | PVC-cable 4-pin, length 5 m, 10 m, 25 m |
| Programming adap | oter |
| MPU-P 9701 | Programming adapter for MPU-4, MPU-H and MPU-M |
| Sealing ring | |

DRE-5Sealing ring for EPA-8, Ø 5 x 1.5 mm, material EPDM (FDA compliant, USP class VI)DRE-15Sealing ring for EPA-18, Ø 15 x 1.5 mm, material EPDM (FDA compliant, USP class VI)









Order code for version with 1x Pt100

| TFP-641 TFP-661 TFP-881 TFP-861 TFP-881 | (for PHARMa (for PHARMa (for PHARMa (for PHARMa (for PHARMa (for PHARMa Sensor Leng 010 025 050 100 Sensor Leng 020 050 | adapt EPA-8, connection head Ø 49 mm, non-sensitive to vibrations) adapt EPA-8, connection head Ø 18 mm, electrical connection via M12 plug) adapt EPA-8, connection head Ø 49 mm, non-sensitive to vibrations) adapt EPA-18, connection head Ø 49 mm, non-sensitive to vibrations) adapt EPA-18, connection head Ø 18 mm, electrical connection via M12 plug) adapt EPA-18, connection head Ø 18 mm, electrical connection via 2.5 m PTFE-cable) sth EL for TFP-6xx in mm (length 10 mm) (length 50 mm) (length 100 mm) sth EL for TFP-8xx in mm (length 100 mm) (length 50 mm) (length 50 mm) Accuracy Class Pt100 A 1/3B | | | | | | | |
|---|--|---|---|--|--|--|--|--|--|
| | | | Electrical PG M12 | Electrical Connection for TFP-641 and TFP-841 PG (cable gland M16x1.5) M12 (M12 plug, standard with MPU-LCD) | | | | | |
| | | | | Transmitter | | | | | |
| | | | | х | (without) | | | | |
| | | | Transmitter for TFP-641 and TFP-841MPU-4(programmable)MPU-H(HART-protocol)MPU-LCD(with display) | | | | | | |
| | | | | Transmitter for TFP-661 and TFP-861 | | | | | |
| | | | | MPU-M | (programmable | 2) | | | |
| | * | | | | Measurement transmitter; no -1040 050 0100 0150 0200 xxyy | Range(only for types with ot selectable at MPU-LCD (measuring range -10+40 °C) (measuring range 0+50 °C) (measuring range 0+150 °C) (measuring range 0+200 °C) (special range) | | | |
| TFP-641/ | 025 / | Α/ | M12 / | MPU-4 / | 0100 | | | | |



The clamp ring is not included in scope of delivery and must be ordered separately.

SRC-05Clamp-tension ring for EPA-8, material 1.4301SRC-10Clamp-tension ring for EPA-18, material 1.4301

Build-In systems

Suitable build-in systems for temperature sensors TFP-641, -642, -661, -681, -841, -842, -861, -881 you will find in product information **Process Connection PHARMadapt EPA**.

| Order code for version with 2x Pt100 | | | | | | | | | | |
|--------------------------------------|---|---|----------------------------|--|-----------------------|--------------|----------------|--------------------------------------|--|--|
| TFP-641.2 | (for PHAR) | Madapt EPA | -8, connectio | on head Ø 4 | 49 mm, 2x P | t100, non-se | ensitive to vi | brations, | | |
| TFP-642 | (for PHARMadapt EPA-8, like TFP-641.2, but higher connection head and prepared for 2x transmitter) | | | | | | | | | |
| TFP-841.2 | (for PHARMadapt EPA-18, connection head Ø 49 mm, 2x Pt100, non-sensitive to vibrations, no transmitter possible!) | | | | | | | | | |
| TFP-842 | (for PHARMadapt EPA-18, like TFP-841.2, but higher connection head and prepared for 2x transmitter) | | | | | | | | | |
| | Sensor Le | ngth EL for | TFP-6xx in | mm | | | | | | |
| | 010 025 | (length 1 | 0 mm) 5 mm) | | | | | | | |
| | 050 | (length 5 | 0 mm) | | | | | | | |
| | Sensor Length EL for TFP-8xx in mm | | | | | | | | | |
| | 020 (length 20 mm) 050 (length 50 mm) | | | | | | | | | |
| | | Accuracy | Class Pt100 |) | | | | | | |
| | | A 1 /7 P | | | | | | | | |
| | | 1/3B 1/10B | | | | | | | | |
| | | | Electrical | Connectio | on for TFP-6 | 41.2 and TF | P-841.2 | | | |
| | | | PG | (cable gl | and M16x1 | .5) | | | | |
| | | | 2xPG 2xM12 | (2x Cable (2x M12 | e gland Mito plug) | X1.5) | | | | |
| | | Electrical Connection for TFP-642 and TFP-842 | | | | | | | | |
| | | M12 (M12 plug) | | | | | | | | |
| | | | 2xM12 (2x M12 plug) | | | | | | | |
| | | Continue if TFP-642 or TFP-842 is selected! No further options for TFP-641.2 or TFP-841.2! | | | | | | | | |
| | | 1. Transmitter | | | | | | | | |
| | | | | MPU-4 (programmable) | | | | | | |
| | | | | Measurement Range 1. MPU | | | | | | |
| | | | | -1040 (measurement range $-10+40$ °C) | | | | | | |
| | | | | 0100 (measurement range 0+100 °C) | | | | | | |
| | | | | 0150 (measurement range $0+150$ °C) | | | | | | |
| | | | | xxyy (special range) | | | | | | |
| | | | | | | 2. Transm | itter | | | |
| | | | | | | MPU-4 | (programr | nable) | | |
| | | | | | | | Measurem | nent Range 2. MPU | | |
| | | | | | | | -1040 050 | (range -10+40 °C) (range 0+50 °C) | | |
| | | | | | | | 0100 | (range 0+100 °C) | | |
| | | | | | | | 0150 0200 | (range 0+150 °C) (range 0+200 °C) | | |
| | | | | | | | ххуу | (special range) | | |
| ♥ TFP-642 / | ¥ 025 / | ¥ A / | ∀ M12 / | ♦ MPU-4.4 | ∜ 050 / | ¥ MPU-4 / | ¥ 050 | | | |
| | 0131 | ~ | 11221 | | 0111907 | | 0 | | | |